TEST DE ANTÍGENOS AUTODIAGNÓSTICO COVID-19

CE
0123

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Product Service

EC Certificate

EC Design-Examination Certificate
Directive 98/79/EC on In Vitro Diagnostic Medical Devices (IVDD), Annex III (6)
(Devices for self-testing)
No. V9 061317 0006 Rev. 00

Manufacturer: **Boson**
90-94 Tianfeng Road
Jimei North Industrial Park
361021 Xiamen, Fujian
PEOPLE'S REPUBLIC OF CHINA

Product: **In Vitro diagnostic devices for self testing**

Model(s): **Rapid SARS-CoV-2 Antigen Test Card**

Parameters:

Model Name:	Model No.:
Rapid SARS-CoV-2 Antigen Test Card	REF 1N40C5-2
Rapid SARS-CoV-2 Antigen Test Card	REF 1N40C5-4
Rapid SARS-CoV-2 Antigen Test Card	REF 1N40C5-6

The Certification Body of TÜV SÜD Product Service GmbH declares that a design examination has been carried out on the respective devices in accordance with IVDD Annex III (6). The design of the devices conforms to the requirements of this Directive. All applicable requirements of the testing and certification regulation of TÜV SÜD Group have to be complied with. For details and certificate validity see: www.tuvsud.com/ps-cert?q=cert:V9_061317_0006_Rev_00

Report No.: 713210321

Valid from: 2021-04-01
Valid until: 2022-05-26

Date, 2021-04-01

Christoph Dicks
Head of Certification/Notified Body

Page 1 of 1
TÜV SÜD Product Service GmbH is Notified Body with identification no. 0123

TÜV SÜD Product Service GmbH • Certification Body • Ridlerstraße 65 • 80339 Munich • Germany

SIN RECETA MÉDICA

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BOSON	
Declaration of Conformity	
Manufacturer	Boson 90-94 Tianfeng Road, Jimei North Industrial Park, Xiamen, Fujian 361021, P. R. China.
European Representative	Lotus NL B.V. Koningin Julianaplein 10, 1e Verd, 2595AA, The Hague, Netherlands.
Product	Rapid SARS-CoV-2 Antigen Test Card, REF: 1N40C5-2 Rapid SARS-CoV-2 Antigen Test Card, REF: 1N40C5-4 Rapid SARS-CoV-2 Antigen Test Card, REF: 1N40C5-6
Classification	Self-Test IVD ←
Confirmatory Assessment Route	Directive 98/79/EC on In Vitro Diagnostic Medical Devices (IVDD), Annex III section 6
We herewith declare that above mentioned products meet the provisions of the council Directive 98/79/EC for medical devices. All Supporting Documentation is Retained under the premises of the manufacturer. We have sole responsibility for issuing the Declaration of Conformity.	
Standard Applied	ISO 13485:2016 EN ISO 14971:2012 EN ISO 18113-1:2011 EN ISO 18113-4:2011 ISO 15223-1:2016 EN 62366-1:2015 EN 13612:2002 EN ISO 23640:2015
Notified Body	TÜV SÜD Product Service GmbH, Ridlerstraße 65, 80339 Munich, Germany. Identification number: 0123
(EC) Certificate(s)	No. V9 061317 0006 Rev. 00
Start of CE-Marking	2021-04-01
Place, Date of Issue	Xiamen, 2021-04-01
Signature	<i>Changgong Zhang</i> (Signed By Boson Representative) Name: Changgong Zhang Title: General Manager
CE 0123	

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