Pedalier Electrónico Loovo

English

Use instruction:

1. Take pedal exerciser out of box, then handle it with care. Look at picture 1.

2. Lay out pedal exerciser on the ground in the same level, then put knob A onto it as picture shows. Look at picture 2.

3. Turn knob A so as to make the axis B pressure properly. (Toward left to put on pressure, otherwise decrease pressure). Turn on the calculator button C and rotate axis B, immediately it is running in use. Look at picture 3. (For further details, please read the user manual for the screen).

4. Pull the pin D toward inside after use and fold it well like picture 5 does.

5. Note:
(1) Inspect it and make sure it is placed steadily and in same level before use.
(2) Don't forget to let the pressure be minimum when start to use, then increase it to reach what it does need.

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

Pedalier Electrónico Loovo

USER MANUAL FOR THE SCREEN:


The different screen instructions are explained below.

- **Switch on/Reset:**
This button is used to switch on the device. To do this, simply press it. If you wish to start a new measurement press the button for 5 seconds.
- **Calories:**
Shows the calories which are being burnt.
- **Time:**
Shows the time during which exercise is being done.
- **Turns:**
Refers to the number of turns being made.
- **Revolutions per minute:**
Counts the turns per minute.
- **Scan:**
Enables the time, calories and revolutions per minute to be seen on the screen in a random manner.

Once that the pedal system is switched on, press the "Switch on/Reset" key to be able to view each of the options listed above.

N.B.:

If not in active use, the screen will automatically switch off in a few minutes.

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