



3F, Building Block 2, No. 3400 Gonghexin Road,
Jing'an District - Shanghai 200436, P.R. CHINA
上海市静安区共和新路3400号2幢3层
Tél. : +86 21 68 55 50 32
Fax : +86 21 68 55 50 33
E-mail : ctcshanghai@ctcgroupe.com

TEST REPORT

Report No.: S201010552_1

14 October 2020

huangcongmin@intco.com

APPLICANT: ANHUI INTCO MEDICAL PRODUCTS CO. LTD
(C41483)

(HAITANG ROAD WEST AND YINHUA ROAD
NORTH) SUIXI WUHU MODERN INDUSTRIAL
PARK, SUIXI TOWN
HUAIBEI ANHUI
CHINA

Date of receipt : 14 Aug. 2020
Testing period : 09 Oct. 2020
: 14 Oct. 2020

Buyer: —

Sample description: Disposable Vinyl Gloves (PVC)

Test(s) requested :—
Service :REGULAR
Brand / Section :—
Season :—
End use :—
Factory name :—
Factory code :—

For CE Marking : Yes
Previous report :—
Product category :—
Product type :—
Test stage :FIRST TEST
Supplier name :—
Exported to :—

1. Conclusion:

	Tests description	Conformity
1	Resistance to degradation by chemicals : EN 374-4 : 2019	None

Pass: requirements met Fail: requirements not met None: no requirement for this test N/A: not applicable

Approved by

Henry YAN
Laboratory Manager

The report is issued by CTC Shanghai under its General Conditions printed overleaf. The results shown in this report refer only to the sample (s) tested. Except by special arrangement, the test items will not be retained by CTC Shanghai for more than 3 months. The test report shall not be reproduced, except in full, without the written approval of the testing laboratory.

To declare the conformity to the requirement, the uncertainty of measurement, associated to the test results, has not been taken into account.



TEST REPORT

Report No. : CH:TX:1142011143

DATE : 19/03/2020




CASE NO: QDHL2002000966MD

ANHUI INTOCO MEDICAL PRODUCTS CO, LTD.
 NO.1, HAITANG ROAD, SUIXI DISTRICT ECONOMIC DEVELOPMENT AREA, HUAIBEI CITY,
 ANHUI PROVINCE.
A/C F619701 SGS-CSTC STANDARD TECHNICAL SERVICES (QINGDAO) CO., LTD.
CONTACT PERSON : --

THE FOLLOWING SAMPLE(S) WAS/WERE SUBMITTED AND IDENTIFIED BY/ON BEHALF OF THE CUSTOMER AS :

SAMPLE DESCRIPTION GLOVES
 DISPOSABLE VINYL GLOVE(PVC)
QD REFERENCE NO QDHL2002000966MD
PHOTO APPENDIX.



SAMPLE RECD ON 03/03/2020 **TESTING PERIOD :** 03/03/2020 - 19/03/2020

Summary of Test Results/Conclusion

Test Method / Standard	Clause/Test Name	Status / Performance Level
EN 374-4:2013	Resistance to Degradation by Chemicals	
	Sodium hydroxide 40%	Refer results.

Per pro SGS India Private Ltd.

K. PACHAIYAPPAN
ASST. MANAGER

Email your Test Report Related Enquiries at Feedback.SLT@sgs.com

TEST REPORT

Report No. : CH:TX:1142011143

DATE : 19/03/2020



CASE NO: QDHL2002000966MD

RESULTS

EN 374-4:2013 Protective Gloves against Chemicals and Micro Organisms – Determination of resistance to degradation by chemicals

Chemical / CAS NO	Exposure Duration	Test Results		Observation
		Percentage change in puncture resistance		
		Glove sample	Result (%)	
Sodium hydroxide 40% 1310-73-2	60±5 minutes	1	11.2%	No change
		2	15.8%	
		3	12.1%	
		Mean	13.0%	
		Standard Deviation	2.477	

***** End of Report*****

TEST REPORT

Report No.: S201010552_1

14 October 2020

APPLICANT: ANHUI INTCO MEDICAL PRODUCTS CO. LTD
(C41483)

2. Sample(s) description assigned by laboratory:

<u>Size</u>	<u>Analyzed product</u>	<u>Description</u>	<u>Sample information</u>
	GLOVE	white vinyl palm--P	



201010552



TEST REPORT

Report No.: S201010552_1

14 October 2020

APPLICANT: ANHUI INTCO MEDICAL PRODUCTS CO. LTD
(C41483)

3. GLOVE/

white vinyl palm—P

	Method	Client Requirement	Unit	Result	Conformity
▲ Resistance to degradation by chemicals	EN 374-4 : 2019				
Challenge chemical tested				Hydrogene peroxide 30%	
Sample degradation by the chemical tested			%	8.5	
Sample degradation by the chemical tested (2)			%	10.3	
Sample degradation by the chemical tested (3)			%	23.5	
Degradation - Mean			%	14.1	
Degradation - Standard deviation			%	8.2	
Lining				No Lining	
Observation				Swelling	

END OF TEST REPORT

▲: The test was carried out by external accredited laboratory under their accreditation scope.