



Hyperbaric Oxygen Therapy

Henry's Law:

“The amount of dissolved gas in a liquid is proportional to its partial pressure above the liquid.”



The state of gas dissolved in liquid at one ATA pressure.



The state of gas dissolved in liquid at two ATA pressure.

Hyperbaric Oxygen Therapy is the medical use of oxygen at an increased level of pressure, (hyper) means increased, (baric) means pressure. Hyperbaric Oxygen Therapy is a treatment which enhances the body's natural healing process by providing an environment which allows the body to absorb much higher amounts of oxygen than possible at normal atmospheric pressure. It increases the oxygen transport capacity in the blood. This noninvasive therapy is the most trusted way to increase oxygen levels to all organs of the body.

Applications of hyperbaric chambers



Adjuvant treatment of some diseases



Rapid recovery after exercise




Family health management



Beauty salon anti-aging

ST1700 Product Information



Medium size sitting type for 1 person, come with a folding chair.

Chamber Parameters

- Size(L*W*H): 170*70*110 cm/ 67*28*44"
- Application: Clinic, Gym, Spa, Homecare
- Pressure: 1.3 ATA/ 6 PSI
- Package dimension: 87*58*47 cm
- Net weight: 33 kg
- Gross weight: 34 kg



Oxygen Concentrator

Imported grade high performance molecular sieve Effectively filter 99% of bacteria Provide you with ultra-pure forest-like oxygen.

- Size(L*W*H): 35*10*65cm / 14*4*26"
- Flow: 1-5L/min (adjustable)
- Power: 500W



Air Cooler

Sufficient power, built-in filter cotton to purify the compressed air, equipped with a silencer to minimise noise.

- Size(L*W*H): 18*12*35cm / 7.2*5*14"
- Power: 200W



Air Compressor


Sufficient power, built-in filter cotton to purify the compressed air, equipped with a silencer to minimise noise.

- Size(L*W*H): 39*24*26cm / 16*10*10"
- Flow: 72L/Min
- Power: 550W

ST1700 Features

- Suitable size for 1 person
- Designed to add a folding chair(adjustable height) for ultimate comfort.
- Can use a laptop or do some reading inside
- Internal metal frame - Hold the shape when deflated.
- Chamber protection cover, to avoid dirty and easy to wash.
- Emergency valve - To speed the depressurization in an emergency.
- Easy Operation - One person can operate it without assistance.
- Delivers 93% oxygen under pressure via an oxygen headset/facial mask.
- Internal Pressure Gauge - User can observe the pressure while inside the chamber.
- Pressure options: 1.3/1.4/1.5ATA

Product Superiority



Cabin raw materials

TPU material passed SGS and Rohs environmental certification, high strength non-toxic environmentally Good elasticity and abrasion resistance, tear resistance, high tensile strength.

Welding

Energy-saving and environmentally friendly high-frequency welding a molding process safe and more stable.



Optional outer cover

Half wrap and full wrap cloth covers are available for choose Different Color can be chosen 100% cotton stretch fabric. High wear-resistant and dust-resistant polyamide fiber cloth cover



Zipper

Yoshida Kogyo Kabushikigaisha The world's top YKK brand smooth zipper The originator of the zipper industry, representing the industry standard Adopt Japanese precise process.



3 types of safety valves ensure safe and effective use



Pressure Gauges



Automatic Pressurize Relief Valves



Emergency Relief Valve



Air Deflate Value

CONTACT US



Whatapp/Tel: +86 18017023916

STS Web: www.macypan.com

Email: Rachel@macy-pan.com



Address: Bldg 15, No.889 Guinan Road, Songjiang District, Shanghai, China 201617

Product Pictures & Customer Feedback

